


LINKLOOP – DEXCOM DATA FLOW DIAGRAM

Partner: VibeCMD LLC | App: LinkLoop | Date: February 2026

ACTORS & SYSTEMS

[USER]	T1D Warrior using LinkLoop iOS app
[APP]	LinkLoop iOS (React Native / Expo, iPhone only)
[SERVER]	LinkLoop Backend (Node.js / Express on Render.com)
[DB]	MongoDB Atlas (cloud database, US region)
[DEXCOM]	Dexcom API – api.dexcom.com (v3 REST API)
[CIRCLE]	Care Circle Members (family / caregivers)

PHASE 1 – DEXCOM AUTHORIZATION (One-time setup, explicit user action)



Saves tokens to [DB] – server-side only, NEVER sent to app

```
db.users.dexcom = {
  accessToken: "...", ← encrypted at rest in MongoDB Atlas
  refreshToken: "...", ← encrypted at rest in MongoDB Atlas
  tokenExpiry: DateTime,
  connected: true,
  lastSync: null
}
```

Returns success HTML page ("Dexcom connected! Close this window.")

[USER] closes browser, returns to LinkLoop app

App polls /api/dexcom/status – sees connected: true

CGM screen now shows "Dexcom Connected" status

PHASE 2 – DATA SYNC (User-initiated, on demand)

[USER]

Taps "Sync Now" on CGM screen

[APP] — HTTPS POST /api/dexcom/sync (JWT auth header) —→ [SERVER]

token expiry:

min remaining → auto-refresh

/v3/oauth2/token (refresh_token)

s new tokens to [DB]

| Checks

| If < 5

| POST

| Save

| Fetches

s glucose readings:

[DEXCOM API] ← HTTPS GET /v3/users/self/egvs —→ [SERVER]

params: startDate, endDate

header: Authorization: Bearer {access_token}

Returns EGV records:

- systemTime (timestamp)
- value (mg/dL integer)
- trend (e.g. "singleUp", "flat")

[SERVER] ← [DEXCOM API]

Transforms & stores in [DB] – GlucoseReading documents:

```
{
  userId: linkloop_user_id,
  value: integer (mg/dL),
  trend: "rising" | "stable" | "falling" | etc,
  trendArrow: "↑" | "→" | "↓" | etc,
  source: "dexcom",
  timestamp: DateTime,
  unit: "mg/dL"
}
```

NOT stored: raw Dexcom userId, personal identifiers,
device serial number, raw trend strings

Returns { synced: N, message: "Synced N readings" }

[APP] ← [SERVER]

Fetches latest readings via GET /api/glucose

CGM screen displays:

- Current glucose value + trend arrow
- 24-hour history graph
- Time since last reading
- High / Low / In Range status

PHASE 3 – CARE CIRCLE SHARING (Optional, warrior controls)

[USER / WARRIOR]

Chooses to invite family / caregiver to their Circle
Sends invite code (generated by server)

[CIRCLE MEMBER] joins via invite code in LinkLoop app

Gets linked to warrior's account as role: "member"
Can view warrior's glucose readings – READ ONLY

[APP – Circle Member] — HTTPS GET /api/glucose/member-view → [SERVER]

Verifies circle membership

Returns warrior's readings from [DB]

(same data warrior sees – no extras)

[APP – Circle Member] ← [SERVER]

Circle member sees: glucose value, trend, history
Circle member CANNOT: modify readings, access Dexcom tokens,
see warrior's Dexcom credentials,
share data further

PHASE 4 – DISCONNECTION / REVOCATION

[USER]

Option A: Taps "Disconnect Dexcom" in LinkLoop app

[APP] — HTTPS POST /api/dexcom/disconnect → [SERVER]

ately nulls all tokens in [DB]:	Immedi
om.accessToken = null	dexc
om.refreshToken = null	dexc
om.connected = false	dexc
ng glucose readings are retained	Existi
s own historical data, not Dexcom's)	(user'
	▼
	[DB] tokens
deleted	

ys		Option B: User revokes access directly at dexcom.com
		→ access_token immediately invalidated by Dexcom
		→ Next sync attempt returns 401
		→ Server marks dexcom.connected = false automatically
ys		Option C: User requests full account deletion
		Email: vibetech@vibecmd.net
		→ All user data including glucose readings deleted within 30 da

=====

DATA STORAGE SUMMARY

=====

What is stored in LinkLoop's database (MongoDB Atlas):

Data	Location	Retention
Dexcom access_token	Server DB only	Until disconnect
Dexcom refresh_token	Server DB only	Until disconnect
Glucose readings (value, trend, timestamp)	Server DB	While account active
User account (name, email)	Server DB	While account active
Care Circle membership	Server DB	While account active

What is NEVER stored:

- Dexcom username or password
- Dexcom device serial numbers
- Raw Dexcom internal user IDs
- Any data not directly needed to display glucose to the user

What is NEVER shared:

- Data is never sold to third parties
- No advertising SDKs in the app
- No analytics platforms with access to health data
- Care Circle members only see data the warrior explicitly shares

=====

SECURITY CONTROLS

=====

- ✓ All traffic – HTTPS / TLS 1.2+ enforced, no HTTP fallback
- ✓ Authentication – JWT tokens, short-lived, validated server-side
- ✓ Dexcom tokens – stored server-side only, never transmitted to app
- ✓ Database – MongoDB Atlas with access controls and encryption at rest
- ✓ Server – Render.com managed hosting (automatic security patches)
- ✓ Auto token refresh – silent, 5 min before expiry, no user interruption
- ✓ Disconnect – full token wipe on user request, immediate effect

CONTACT

Company: VibeCMD LLC
App: LinkLoop
Owner: Kevin Cunningham
Email: vibetech@vibecmd.net
Privacy: <https://vibecmd.net/privacy/linkloop>
Server: <https://linkloop-9l3x.onrender.com>
